

## Anti-Gastrin 17 antibody (STJ401381)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/CLIA/ELISA/IFA
Host / Source	Mouse
Reactivity	G-17

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein A Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Storage Instruction	Short Term ( $\leq 2$ weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	Gastrin 17.
-------------	-------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	-------------------------------------------------------------------------------------------------------------------------