

Anti-Feline Immunodeficiency Virus antibody (STJ401385)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ELISA/IFA
Host / Source	Mouse
Reactivity	Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.2
Isotype	IgG1
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	FIV p24 (mature protein). No reactivity with FeLV, Feline Herpesvirus type 1, Feline Coronavirus, Feline Calicivirus
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---