

## Anti-Progastrin-releasing Peptide antibody (STJ401389)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | CLIA/ELISA/CMA     |
| Host / Source | Mouse              |
| Reactivity    | ProGRP             |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Monoclonal   |
| Conjugation         | Unconjugated   |
| Purification        | Purified Monoclonal Antibody   |
| Formulation         | 10 mM Phosphate Buffer with 150 mM NaCl, pH 7.4  |
| Isotype             | IgG2a  |
| Storage Instruction | Short Term ( $\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles. |

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|