

Anti-Campylobacter antibody (STJ401404)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	Campylobacter

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	Phosphate Buffered Saline, pH 7.4
Isotype	IgG1
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	No reactivity with <i>H. pylori</i> , <i>E. coli</i> , <i>S. enteritidis</i> , <i>S. paratyphi</i> , <i>S. typhi</i> , <i>S. typhimurium</i> , <i>S. boydii</i> , <i>S. dysenteriae</i> , <i>S. flexneri</i> , <i>S. sonnei</i> , <i>Y. enterocolitica</i> , <i>L. monocytogenes</i> , <i>S. aureus</i>
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---