

Anti-Haptoglobin antibody (STJ401431)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/LF/ELISA
Host / Source	Mouse
Reactivity	Cat

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.2
Storage Instruction	Store at 2–8 °C, and do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---