

Anti-6-Monoacetylmorphine antibody (STJ401433)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	6-MAM

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Isotype	IgG1k
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	6-Monoacetylmorphine and Diacetylmorphine. Morphine, Codeine, Hydrocodone, Hydromorphone, Oxycodone, Ethylmorphine, Oxymorphone, and Hydromorphone = 1%
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---