

Anti-Bovine Growth Hormone antibody (STJ401435)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Mouse
Reactivity	rBGH

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Ion Exchange Purified Monoclonal Antibody
Formulation	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Isotype	IgG2a
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Recombinant Bovine Growth Hormone. Pituitary Bovine Growth Hormone (pBGH) , Native Bovine Growth Hormone (nBGH)
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---