

Anti-CYFRA21-1 antibody (STJ401469)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CLIA/ELISA
Host / Source	Mouse
Reactivity	Cytokeratin Fragment

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Isotype	IgG2a
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---