

Anti-Mouse Immunoglobulin antibody (STJ401505)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Goat
Reactivity	IgG

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Immunoaffinity purified Polyclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Storage Instruction	Store at 2–8 °C, and do not freeze.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---