

Anti-Bovine Serum Albumin antibody (STJ401507)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC/ELISA
Host / Source	Rabbit
Reactivity	BSA

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified Polyclonal Antibody
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Single precipitin arc using immunoelectrophoresis against BSA and the albumin in bovine serum.
-------------	--

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---