

Anti-Apolipoprotein A1 antibody (STJ401513)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	TIA
Host / Source	Goat
Reactivity	ApoA1

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Polyclonal, Delipidated and 0.2 μ m filtered
Formulation	20 mM Tris at pH 7.4 with 250mM NaCl
Storage Instruction	Short Term ($\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Single precipitin arc visible with immunoelectrophoresis.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---