

## Anti-Apolipoprotein Lp (a) antibody (STJ401515)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | TIA                |
| Host / Source | Goat               |
| Reactivity    | LpA                |

### PRODUCT PROPERTIES

|                     |                                                                                              |
|---------------------|----------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                   |
| Conjugation         | Unconjugated                                                                                 |
| Purification        | Polyclonal, Delipidated and 0.2 $\mu$ m filtered                                             |
| Formulation         | 0.05M Tris, pH 7.5 with 0.5M NaCl                                                            |
| Storage Instruction | Short Term ( $\leq$ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |                                                           |
|-------------|-----------------------------------------------------------|
| Specificity | Single precipitin arc visible with immunoelectrophoresis. |
|-------------|-----------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|