

Anti-Ceruloplasmin antibody (STJ401519)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IEP/TIA
Host / Source	Goat
Reactivity	Ceruloplasmin

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Polyclonal Antibody, IgG Fraction
Formulation	20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Human Ceruloplasmin.
-------------	----------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---