

Anti-6-Monoacetylmorphine antibody (STJ401521)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Rabbit
Reactivity	6-MAM

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Protein A Purified Polyclonal Antibody, Ig Fraction
Formulation	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	6-Monoacetylmorphine. Desomorphine
-------------	------------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---