

## Anti-IgE antibody (STJ401523)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/ELISA       |
| Host / Source | Goat               |
| Reactivity    | Human              |

### PRODUCT PROPERTIES

|                     |                                                            |
|---------------------|------------------------------------------------------------|
| Clonality           | Polyclonal                                                 |
| Conjugation         | Unconjugated                                               |
| Purification        | Affinity Purified Polyclonal Antibody                      |
| Formulation         | 10 mM Sodium Phosphate, pH 7.2 with 0.15 M Sodium Chloride |
| Storage Instruction | Store at 2–8 °C, and do not freeze.                        |

### TARGET INFORMATION

|             |                                                                                                               |
|-------------|---------------------------------------------------------------------------------------------------------------|
| Specificity | Human IgE Heavy (e) chain. No reactivity with Human immunoglobulin light chains or other Human serum proteins |
|-------------|---------------------------------------------------------------------------------------------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|