

Anti-IgA antibody (STJ401530)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/TIA
Host / Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Polyclonal, Delipidated and 0.2 μ m filtered
Formulation	20 mM Sodium Phosphate, pH 7.2 with 150 mM Sodium Chloride
Isotype	IgG
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	Single precipitin arc visible with immunoelectrophoresis.
-------------	---

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---