

Anti-Histone Core antibody (STJ401532)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Sheep
Reactivity	Histone Core

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified Polyclonal Antibody
Formulation	10 mM Phosphate Buffered Saline
Storage Instruction	Store at-20 °C and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Reactive with many species
-------------	----------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---