

Anti-Histone H1 antibody (STJ401533)

GENERAL INFORMATION

| | |
|---------------|--------------------|
| Product Type | Primary antibodies |
| Applications | WB/IHC/IP |
| Host / Source | Sheep |
| Reactivity | Histone H1 |

PRODUCT PROPERTIES

| | |
|---------------------|---|
| Clonality | Polyclonal |
| Conjugation | Unconjugated |
| Purification | Purified Polyclonal Antibody |
| Formulation | 10 mM Phosphate Buffered Saline |
| Storage Instruction | Store at-20 °C and avoid repeated freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------|----------------------------|
| Specificity | Reactive with many species |
|-------------|----------------------------|

ADDITIONAL INFORMATION

| | |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|