

Anti-Histone H2B antibody (STJ401534)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IP
Host / Source	Sheep
Reactivity	Histone H2B

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Purified Polyclonal Antibody
Formulation	10 mM Phosphate Buffered Saline
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Reactive with many species
-------------	----------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---