

## Anti-Histone H3 antibody (STJ401535)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IHC/IP          |
| Host / Source | Sheep              |
| Reactivity    | Histone H3         |

### PRODUCT PROPERTIES

|                     |                                                        |
|---------------------|--------------------------------------------------------|
| Clonality           | Polyclonal                                             |
| Conjugation         | Unconjugated                                           |
| Purification        | Purified Polyclonal Antibody                           |
| Formulation         | 10 mM Phosphate Buffered Saline                        |
| Storage Instruction | Store at -20 °C and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |                            |
|-------------|----------------------------|
| Specificity | Reactive with many species |
|-------------|----------------------------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|