

Anti-Prealbumin antibody (STJ401543)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Goat
Reactivity	Prealbumin

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Polyclonal Antibody, defibrinated, delipidated, absorbed by solid phase chromatography as required and 0.2 μ m filtered
Formulation	0.05M Tris, pH 7.5 with 0.5M NaCl
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Human.
-------------	--------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---