

## Anti-Apolipoprotein CIII antibody (STJ401547)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Goat
Reactivity	Apo CIII

### PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity Purified Polyclonal Antibody
Formulation	75 mM Sodium Phosphate, 75 mM NaCl, 0.5 mM EDTA, pH 7.2
Storage Instruction	Short Term ( $\leq 2$ weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	Human Apo CIII.
-------------	-----------------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---