

Donkey Anti-Sheep IgG H&L antibody {Unconjugated} (STJ401549)

GENERAL INFORMATION

Product Type	Secondary antibodies
Applications	ELISA
Host / Source	Donkey
Reactivity	Sheep

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Affinity Purified Polyclonal Antibody
Formulation	PBS, pH 7.2
Isotype	IgG
Storage Instruction	Store at 2–8 °C, and do not freeze.

TARGET INFORMATION

Specificity	H & L chains Sheep IgG.
-------------	-------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---