

Anti-Carisoprodol antibody (STJ401558)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Rabbit
Reactivity	Soma

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Protein A Purified Polyclonal Antibody, Ig Fraction
Formulation	0.05 M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	Carisoprodol. meprobamate 12%
-------------	-------------------------------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---