

## Anti-Herpes Simplex Virus 1 antibody (STJ401573)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IFA
Host / Source	Mouse
Reactivity	HSV 1-gE

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG
Storage Instruction	Store at -20 °C and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	gE-1.
-------------	-------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---