

## Anti-Herpes Simplex Virus 2 antibody (STJ401576)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IFA             |
| Host / Source | Mouse              |
| Reactivity    | HSV 2-ICP27        |

### PRODUCT PROPERTIES

|                     |                                                        |
|---------------------|--------------------------------------------------------|
| Clonality           | Monoclonal                                             |
| Conjugation         | Unconjugated                                           |
| Purification        | Protein G Purified Monoclonal Antibody                 |
| Formulation         | 10 mM Phosphate Buffered Saline, pH 7.4                |
| Isotype             | IgG                                                    |
| Storage Instruction | Store at -20 °C and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|             |        |
|-------------|--------|
| Specificity | ICP27. |
|-------------|--------|

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|