

Anti-Herpes Simplex Virus 2 antibody (STJ401580)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IFA
Host / Source	Mouse
Reactivity	HSV 2-gC

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG2ak
Storage Instruction	Store at-20 °C and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	gC-2.
-------------	-------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---