

## Anti-Cytomegalovirus antibody (STJ401585)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IFA
Host / Source	Mouse
Reactivity	Cytomegalovirus (CMV-pp28)

### PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG2ak
Storage Instruction	Store at -20 °C and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

Specificity	pp28.
-------------	-------

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---