

Anti-Cytomegalovirus antibody (STJ401586)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IFA
Host / Source	Mouse
Reactivity	CMV-IE1/2

PRODUCT PROPERTIES

Clonality	Monoclonal
Conjugation	Unconjugated
Purification	Protein G Purified Monoclonal Antibody
Formulation	10 mM Phosphate Buffered Saline, pH 7.4
Isotype	IgG1k
Storage Instruction	Store at -20 °C and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	IE1/2.
-------------	--------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---