

Anti-Cannabinoid antibody (STJ401589)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Sheep
Reactivity	THC

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	DEAE Purified Polyclonal Antibody; Ig Fraction
Formulation	20 mM Phosphate Buffer, pH 7.1 with 150 mM Sodium Chloride
Storage Instruction	Short Term (\leq 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Specificity	THC.
-------------	------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---