

Anti-Benzodiazepine antibody (STJ401591)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Sheep
Reactivity	Benzodiazepine (Oxazepam)

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	DEAE Purified Polyclonal Antibody, Ig Fraction
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---