

Anti-Cotinine-4 antibody (STJ401597)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	LF/ELISA
Host / Source	Rabbit
Reactivity	Cotinine-4

PRODUCT PROPERTIES

Clonality	Polyclonal
Conjugation	Unconjugated
Purification	Protein A Purified Polyclonal Antibody, Ig Fraction
Formulation	10 mM Phosphate Buffered Saline, pH 7.0
Storage Instruction	Short Term (≤ 2 weeks) : +2–8 °C. Long Term:-20 °C. Avoid repeated freeze-thaw cycles.

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---