

Anti-ACTH antibody (138-176) {FITC} (STJ500038)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CM/DB/ELISA/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.63-0.68 $\mu\text{g}/\mu\text{l}$
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:1, 000 IP: 1:200 IHC: 1:250 ICC: 1:250 IF: 1:200 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	138-176
------------------	---------

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Slot blot with Anti-ACTH Antibody (FITC) (STJ500038) of rat at 1:500 dilution in buffer. Varying concentration of ACTH was blotted on Nitrocellulose.