

## Anti-AMACR antibody (260-310) (STJ500070)

### GENERAL INFORMATION

|               |                              |
|---------------|------------------------------|
| Product Type  | Primary antibodies           |
| Applications  | CM/DB/ELISA/ICC/IF/IHC/IP/WB |
| Host / Source | Rabbit                       |
| Reactivity    | Human/Mouse                  |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                           |
| Concentration       | 0.6 µg/µl                                                                                                                                                            |
| Conjugation         | Unconjugated                                                                                                                                                         |
| Purification        | Affinity Purified                                                                                                                                                    |
| Dilution Range      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:200<br>IHC: 1:100<br>ICC: 1:100<br>IF: 1:200<br>CFM: 1:100                                                   |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG                                                                                                                                                                  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">23600</a>       |
| Gene Symbol      | <a href="#">AMACR</a>       |
| UniProt ID       | <a href="#">AMACR_HUMAN</a> |
| Immunogen Region | 260-310                     |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot with Anti-AMACR Antibody (STJ500070) of HEK293 cell at 1:500 dilution in buffer.