

Anti-AMACR antibody (260-310) {FITC} (STJ500072)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	CM/DB/ELISA/ICC/IF/IHC/IP/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.6 µg/µl
	Conjugation	FITC
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200 IHC: 1:100 ICC: 1:100 IF: 1:200 CFM: 1:100
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	23600
	Gene Symbol	AMACR
	UniProt ID	AMACR_HUMAN
	Immunogen Region	260-310
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.