

Anti-APH1B antibody (235-257) {Biotin} (STJ500106)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	C.elegans/Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	0.6 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	83464
Gene Symbol	APH1B
UniProt ID	APH1B_HUMAN
Immunogen Region	235-257
Specificity	Does not cross react to other components of gamma secretase complex nor with other proteins.
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.