

Anti-SLC7A10 antibody (1-50) {Biotin} (STJ500169)

GENERAL INFORMATION

| | |
|---------------|---------------------------|
| Product Type | Primary antibodies |
| Applications | CM/ELISA/ICC/IF/IHC/IP/WB |
| Host / Source | Rabbit |
| Reactivity | Mouse/Rat |

PRODUCT PROPERTIES

| | |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality | Polyclonal |
| Concentration | 0.50-0.75 $\mu\text{g}/\mu\text{l}$ |
| Conjugation | Biotin |
| Purification | Affinity Purified |
| Dilution Range | WB: 1:500-1:1,000 DB: 1:4,000 ELISA: 1:4,000 IP: 1:200 IHC: 1:50 ICC: 1:50 IF: 1:50 CFM: 1:50 |
| Formulation | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide). |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|------------------|----------------------------|
| Gene ID | 53896 |
| Gene Symbol | Slc7a10 |
| UniProt ID | AAA1_MOUSE |
| Immunogen Region | 1-50 |

ADDITIONAL INFORMATION

| | |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|