

## Anti-BEST4 antibody (400-500) (STJ500238)

### GENERAL INFORMATION

|               |                      |
|---------------|----------------------|
| Product Type  | Primary antibodies   |
| Applications  | ELISA/IP/WB          |
| Host / Source | Rabbit               |
| Reactivity    | D.melanogaster/Human |

### PRODUCT PROPERTIES

|                     |                                                                                                                                                                      |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                           |
| Concentration       | 0.75-1.20 µg/µl                                                                                                                                                      |
| Conjugation         | Unconjugated                                                                                                                                                         |
| Purification        | Affinity Purified                                                                                                                                                    |
| Dilution Range      | WB: 1:500-1:1, 000<br>DB: 1:2, 000<br>ELISA: 1:2, 000<br>IP: 1:200                                                                                                   |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG                                                                                                                                                                  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.                                                                                                      |

### TARGET INFORMATION

|                  |                             |
|------------------|-----------------------------|
| Gene ID          | <a href="#">266675</a>      |
| Gene Symbol      | <a href="#">BEST4</a>       |
| UniProt ID       | <a href="#">BEST4_HUMAN</a> |
| Immunogen Region | 400-500                     |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|