

Anti-Calpain antibody (700-758) (STJ500354)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Rabbit
Reactivity	S.mansoni

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.65 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	700-758
------------------	---------

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot with Anti-Calpain Antibody (STJ500354) of rat at 1: dilution in buffer. Purified rsm-p80 (concentration 0.31 mg/mL) from E. Coli, 3µg sample was loaded on SDS-PAGE and then transferred to nitrocellulose membrane. Antibodies: Rabbit Anti-sm-p80 Antibody (Calpain), HRP anti-rabbit IgG 1:2000