

## Anti-CCR4 antibody (1-50) {FITC} (STJ500407)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	CF/DB/ELISA/IF/IHC/IP/WB
Host / Source	Rabbit
Reactivity	Cat/Human/Monkey/Mouse/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.68 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 IP: 1:250 IHC: 1:250 ICC: 1:250 IF: 1:250 CFM: 1:250
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">1233</a>
Gene Symbol	<a href="#">CCR4</a>
UniProt ID	<a href="#">CCR4_HUMAN</a>
Immunogen Region	1-50

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---