

Anti-CHN2 antibody (1-100) {Biotin} (STJ500532)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/WB
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.68-0.71 $\mu\text{g}/\mu\text{l}$
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000 IP: 0 IHC: 0 ICC: 0 CFM: 0
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at-20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	1124
	Gene Symbol	CHN2
	UniProt ID	CHIO_HUMAN
	Immunogen Region	1-100
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.