

**Anti-CCR2 antibody (C-Term) {FITC} (STJ500653)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	CM/ELISA/ICC/IF/IHC/IP/FC/WB
Host / Source	Rabbit
Reactivity	Human

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Concentration	0.50-0.60 $\mu\text{g}/\mu\text{l}$
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10,000 ELISA: 1:10,000 IP: 1:250 IHC: 1:200 ICC: 1:200 IF: 1:200 CFM: 1:200
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">729230</a>
Gene Symbol	<a href="#">CCR2</a>
UniProt ID	<a href="#">CCR2_HUMAN</a>
Immunogen Region	C-Term

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---