

## Anti-DPYSL4 antibody (509-527) {FITC} (STJ500797)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Cat/Human/Monkey/Mouse/Pig/Rat

### PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.55 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:450 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">10570</a>
Gene Symbol	<a href="#">DPYSL4</a>
UniProt ID	<a href="#">DPYL4_HUMAN</a>
Immunogen Region	509-527

### ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---