

Anti-ELOVL4 antibody (294-301) {Biotin} (STJ500877)**GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Rabbit
Reactivity	Human/Monkey/Pig/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.6 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: 1:500
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6785
Gene Symbol	ELOVL4
UniProt ID	ELOV4_HUMAN
Immunogen Region	294-301

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---