

## Anti-ergo-1 antibody (790-840) {Biotin} (STJ500958)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | ELISA/IP/WB        |
| Host / Source | Rabbit             |
| Reactivity    | C.elegans          |

### PRODUCT PROPERTIES

|                     |  |
|---------------------|--|
| Clonality           | Polyclonal   |
| Concentration       | 0.62-0.68 $\mu\text{g}/\mu\text{l}$  |
| Conjugation         | Biotin   |
| Purification        | Affinity Purified  |
| Dilution Range      | WB: 1:500<br>DB: 1:10, 000<br>ELISA: 1:10, 000<br>IP: 1:150  |
| Formulation         | Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide). |
| Isotype             | IgG  |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles.  |

### TARGET INFORMATION

|                  |         |
|------------------|---------|
| Immunogen Region | 790-840 |
|------------------|---------|

### ADDITIONAL INFORMATION

|      |   |
|------|---|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|---|