

Anti-FeSODB1 antibody (STJ501002)

GENERAL INFORMATION


Product Type	Primary antibodies
Applications	ELISA/WB
Host / Source	Chicken
Reactivity	L.infantum

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.55 µg/µl
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot with Anti-FeSODB1 Antibody (STJ501002) with recombinant protein. Antibody diluted at 1:500 in buffer. 500 ng and 1 µg protein/lane.