

Anti-C11orf1 antibody (50-100) {Biotin} (STJ501065)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Concentration	0.65-0.67 µg/µl
Conjugation	Biotin
Purification	Affinity Purified
Dilution Range	WB: DB: 1:15, 000-1:20, 000 ELISA: 1:15, 000-1:20, 000
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	64776
Gene Symbol	CFAP68
UniProt ID	CFA68_HUMAN
Immunogen Region	50-100

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---