

Anti-HIPK4 antibody (470-550) (STJ501349)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Rabbit
Reactivity	Human/Monkey/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	0.5-1.5 $\mu\text{g}/\mu\text{l}$
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB: 1:2, 500-1:5, 000 ELISA: 1:50, 000 IP: 1:200 IHC: 1:50-1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgG
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	147746
Gene Symbol	HIPK4
UniProt ID	HIPK4_HUMAN
Immunogen Region	470-550

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Anti-HIPK4 Antibody (STJ501349) of Human at 1:500 dilution in buffer. Lane A) Human HIPK4 over expressing NIH 3T3 cells reduced, Lane B) Mouse HIPK4 over expressing NIH 3T3 cells, Lane C) Mock transfected NIH 3T3 cells. Total protein per lane 2.5 μg .