

**Anti-MAPK3 antibody (300-350) {FITC} (STJ501704)****GENERAL INFORMATION**

Product Type	Primary antibodies
Applications	ELISA/IP/WB
Host / Source	Chicken
Reactivity	Human/Mouse/Rat

**PRODUCT PROPERTIES**

Clonality	Polyclonal
Concentration	0.5 µg/µl
Conjugation	FITC
Purification	Affinity Purified
Dilution Range	WB: 1:500 ELISA: 1:10, 000 IP: 1:150
Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, HEPES, and long-term preservatives (0.02% Sodium Azide).
Isotype	IgY
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

Gene ID	<a href="#">5595</a>
Gene Symbol	<a href="#">MAPK3</a>
UniProt ID	<a href="#">MK03_HUMAN</a>
Immunogen Region	300-350

**ADDITIONAL INFORMATION**

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---