

Anti-MSH4 antibody (1-100) {Biotin} (STJ501806)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	ELISA/WB
	Host / Source	Rabbit
	Reactivity	A.thaliana
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	0.95-1.04 $\mu\text{g}/\mu\text{l}$
	Conjugation	Biotin
	Purification	Affinity Purified
	Dilution Range	WB: 1:500 DB: 1:10, 000 ELISA: 1:10, 000
	Formulation	Contains Tris, HCl/Glycine buffer pH 7.4-7.8, 30% Glycerol and 0.5% BSA, along with cryo-protective agents, Hepes, and long-term preservatives (0.02% Sodium Azide).
	Isotype	IgG
	Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	1-100
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.